

# Blues



Folding Armchair

# Blues – elegance, simplicity, functionality

Blues armchairs not only make sitting comfortable, but can also be used as an extra bed in patients' rooms in hospitals, hospices and so on.



■ Removable cushions increase user comfort.


■ Contoured armrests enhance aesthetic features of the armchair.

■ Luxury version of the armchair comes in genuine leather.

- ELEGANT DESIGN
- MAXIMUM COMFORT
- QUALITY MATERIALS AND DESIGN
- MULTI-PURPOSE USE

## ELEGANCE AND COMFORT

- Blues is an elegant upholstered armchair with a smooth design and folding feature.
- The armchair upholstery may be made of fabric, two types of artificial leather or genuine leather.
- The seat and backrest are made of quality materials for maximum comfort.
- The standard armchair comes with two cushions.
- Once unfolded the head end can be adjusted to two positions.



*Thanks to its colour solution the armchair can become a distinctive design feature in an interior.*

## FUNCTIONS

- The Blues armchair unfolds and may serve as an extra bed.
- The simple unfolding mechanism made of the metal construction ensures stability when lying on the armchair.
- Two castors in the rear corners make the armchair easy to handle.
- The stiffening of the bottom is made of strong beech wood.
- The contoured armrests complete the spacious seating area.



*By unfolding the armchair it is very easy to create a comfortable couch.*

## TECHNICAL PARAMETERS

seat height	50 cm
width	90 cm
height	89 cm
depth	87 cm
unfolded length	203 cm
seating area	67×53 cm
mattress platform (unfolded)	199×67 cm
mattress platform height	41 cm