

## Invitation to EWMA 2016

We would like to invite you to the **European Wound Management Association Conference 2016**, which will be held from the **11th to 13th May 2016 in Bremen, Germany**. At the expert event for healthcare professionals focussing on various aspects of **treatment for wounds and pressure ulcers**, the LINET Group company will present a wide range of innovative solutions aimed at the prevention and treatment of pressure ulcers.



**It will be our great pleasure to welcome you at stand 5E10**, where you can find out about various technologies intended for intensive and long-term care. **We look forward to meeting you!**

### Two professionals: Multicare and Symbioso

The **Symbioso mattress integrated in the Multicare resuscitation bed** is exceptionally silent and effective. The mattress works on the **dual principle of maintaining a constant pressure and suitable microclimate**. The combination of both systems provides effective prevention against the formation of pressure ulcers by eliminating pressure and dampness between the patient's body and the mattress. A suitable microclimatic effect is ensured by the constant airflow in the mattress. In combination with the **positioning options of the Multicare bed**, it has a significant preventative and therapeutic effect.

### Virtuoso: Zero pressure is the best way

You can see the Eleganza 3XC intensive care bed in combination with the Virtuoso active anti-decubitus system. This mattress works on the principle of **alternating 3 linked cells**, which provide **relief and reduction in pressure to zero on the patient's entire body every 7.5 minutes**. The Virtuoso mattress delivers an effective **24-hour care** solution for patients with the greatest risk of pressure ulcer formation.

### Comfortable and safe area

The modern Image 3 bed can be applied universally in hospitals and with a mattress Clinicare 100HF **creates a comfortable and safe area on the bed**. Reactive mattress Clinicare 100HF combines the **benefits of passive and active systems** and offers effective pressure ulcer prevention to a broad spectrum of hospital care. The foam layer creates better pressure relief as opposed to standard mattresses and **the air cells further optimise pressure**.



**EWMA 2016**

11–13 May, 2016

Messe Bremen & ÖVB Arena

Hall 5

Stand no: 5E10

**LINET** Group

LINET Group  
Želečnice 5, 274 01 Slaný, Czech Republic