

ESSENZA ACCESSORIES

LINET




Complete & Seamless Solution

LIFTING POLES & ACCESSORIES


The lifting pole is an accessory used to increase patient independence and mobility.



DESCRIPTION	SWL	COATING	MODEL
Lifting pole to be used with fixed head-end		lacquered RAL 9002	11011410B0000
Lifting pole to be used with fixed head-end	75 kg	lacquered RAL 9006	11011410A0000
Lifting pole to be used with fixed head-end		Chrome plated	1101141000000



DESCRIPTION	SWL	MATERIAL	MODEL
Grey plastic holder to be placed on a lifting pole	75 kg	ABS plastic	4ROTGERSG700-2



DESCRIPTION	SWL	MATERIAL	MODEL
Grey plastic holder to be placed on a lifting pole adjustable by retractor	75 kg	ABS plastic	4ROTGNR100GR-2



DESCRIPTION	SWL	COATING	MODEL
IV holder, 3 hooks to be placed on a lifting pole	9 kg	Powder coated RAL 9002	16010700B0000
		Powder coated RAL 9006	16010700A0000



DESCRIPTION	SWL	COATING	MODEL
IV holder, 3 hooks to be placed on a lifting pole	9 kg	Chrome plated	4MAPL00N1001

IV POLES & ACCESSORIES

Telescopic IV poles can be inserted into the sockets or attached on the bed frame at the head end.



DESCRIPTION	SWL	COATING	MODEL
Telescopic IV pole s-shape, washable	8 kg	Stainless steel	4PV50-0742-0000
Telescopic IV pole s-shape	8 kg	Stainless steel	4MA45MNA0003
Telescopic IV pole s-shape	15 kg	Chrome plated	4PV17604800L



DESCRIPTION	SWL	COATING	MODEL
Telescopic IV pole	15 kg	Chrome plated	4PV17744500L
		Chrome plated	4PV17587200L
		Chrome plated	4PV17744500L
		Chrome plated	4PV297475000
		Stainless steel	4PV322717000
		Chrome plated	4PV340290000



DESCRIPTION	SWL	COATING	MODEL
Holder of IV pumps	3 kg	Chrome plated	1102480000000
			1102520000000



DESCRIPTION	SWL	MATERIAL	MODEL
Infusion bottle basket to be hang on IV pole	Defined by SWL of IV pole	Stainless steel	4DR426101

NAME & CHART HOLDERS

Streamlining patient care with chart and bed holders.



DESCRIPTION	SWL	MATERIAL	MODEL
Name holder for plastic bed end			4GT006000000
Name holder for HPL bed ends	—	Polyethylen	4GT008000000
Name holder for tubular metal or plastic bed ends			4ZZ105313



DESCRIPTION	SWL	COATING / MATERIAL	MODEL
Chart holder for HPL bed ends		RAL 9006	11024001A0000
Chart holder for plastic bed ends	—	ABS plast	11023900A0000

LATERAL TILT SUPPORT

Recommended accessories for making lateralization even safer.



DESCRIPTION	SWL	MATERIAL	MODEL
Positioning cushion for lateral tilt , washable stretch cover (1 piece)	—	Polystyrene beads	4SET17500001



DESCRIPTION	SWL	COATING	MODEL
Ventilation circuit holder	—	Stainless steel	Coming soon

HIGHER SAFETY & COMFORT


Enhance further safety and comfort of the patient.



DESCRIPTION	SWL	COATING / MATERIAL	MODEL
Pair of extenders for collapsible siderails		RAL 9006	11008900A0000
		RAL 9002	11008900B0000
		Aluminium	1100910000000



DESCRIPTION	SWL	MATERIAL	MODEL
Handset holder horizontal for telescopic siderails	—	ABS plastic	1102860000000
Handset holder vertical for telescopic siderails			1102850000000



DESCRIPTION	SWL	MATERIAL	MODEL
Handy box	0.5 kg	Polycarbonate	41029700L0000



DESCRIPTION	SWL	COATING / MATERIAL	MODEL
Crutch holder for plastic and metal bed ends	3 kg	ABS plastic	1102070000000
		RAL 9006	11021100A0000



DESCRIPTION	SWL	MATERIAL	MODEL
Hook for accessories	2 kg	ABS plastic	1102840000000



DESCRIPTION	SWL	COATING	MODEL
Towel holder	3 kg	RAL 9002	11025002B0000
		RAL 9006	11025001A0000



DESCRIPTION	SWL	MATERIAL	MODEL
Plastic siderails covers	—	PUR foam	4PPCHRPEZ300P



DESCRIPTION	SWL	MATERIAL	MODEL
Telescopic siderails covers	—	PUR foam	4PPCHRPEZ300



DESCRIPTION	SWL	MATERIAL	MODEL
Reading lamp	—	Aluminium	4PV383435X00

SHELVES

Shelves and tray to be placed at the bed end for further medical equipment.



DESCRIPTION	SWL	COATING	MODEL
Writing shelf for plastic or metal bed ends	15 kg	RAL 9002	11026500B0000



DESCRIPTION	SWL	MATERIAL	MODEL
Universal utility shelf	10 kg	ABS plastic + iron	11026000B0000

OXYGEN BOTTLE HOLDERS

Oxygen bottle holders can also be fitted to the bed frame.



DESCRIPTION	SWL	MATERIAL	MODEL
Horizontal oxygen bottle (5l) holder to be placed on crossbar on the head end of the bed (Size of O ₂ bottle: Max: L520 × D140; Min: L290 × D90)	15 kg	Stainless steel	4MAR2010PC004



DESCRIPTION	SWL	MATERIAL	MODEL
Vertical oxygen bottle holder to be placed on head end (Size of O ₂ bottle: Max: L800 × D140; Min: L325 × D120)	15 kg	Stainless steel	4MAR250140*
			4MAR250141**

* max bottle diameter 14 cm

** max bottle diameter 12 cm | SWL = Safe Working Load

URINE/DRAINAGE BOTTLES HOLDERS

The stainless steel rail and undercarriage provide hooks and sockets that can be used for urine / drainage bags and bottles, attaching cannulas and more.



DESCRIPTION	SWL	COATING / MATERIAL	MODEL
Urine bottle holder	5 kg	RAL 9002	11027900B0000
Urine sack holder	5 kg	RAL 9002	11023700B0000
Urine bottle metal holder	2 kg	RAL 9006	11027200A0000
Urine bottle holder for telescopic siderails	5 kg	Stainless steel	4MAS6018966



DESCRIPTION	SWL	COATING	MODEL
Redon drainage bottles holder	5 kg	RAL 9002	11008200B0000
Redon drainage bottles holder		RAL 9006	11008200A0000

Please check the availability with your area sales manager or in the configurator.



LINET

Members of LINET Group

LINET spol. s r.o. | Želevčice 5, 274 01 Slaný | Czech Republic
tel.: +420 312 576 400 | fax: +420 312 522 668 | e-mail: info@linet.com | www.linet.com

