


SPRINT 200 & 200 E


LINET

The Stretcher Redefined



Premium Transport & Emergency Stretcher

ERGONOMIC DRIVE THAT GOES THE DISTANCE



i-Drive Power® reduces transport related MSD costs by up to **80%** and drives up to **7 miles** on a single charge.¹



FlexiDrive® makes steering easier for the caregiver and reduces initial push forces by **77%**.¹



IV&Drive® allows the stretcher to be moved from any height without compromising ergonomic posture.¹

SAFETY PROVEN IN MILLION CAREBEDS BROUGHT TO A STRETCHER



Multizone Bed Exit Alarm notifies staff when a patient leaves the stretcher.



Ergoframe® reduces the peak pressure in pelvic area by **16%**.²



Mobi-Lift® provides support for patients as they get off the stretcher.



Scales provide accuracy up to 0.5 kg and eliminate the need of transferring the patient to another surface.

EMPOWERED PATIENT COMFORT, UTMOST EFFICIENCY IN CARE



Electric Positioning reduces physical strain on caregivers during patient repositioning.







Accessories Compatibility ensures streamlined workflows across departments and minimizes issues from incompatibility.






USB Integration allows easy and safe charging of patient's phone.

SPRINT 200 E SPRINT 200

 Support Patient Recovery	<div><div></div></div>	<div><div></div></div>
4-section Mattress Platform	●	●
 Pressure Injury Prevention	<div><div></div></div>	<div><div></div></div>
Ergoframe®	●	●
Advanced Mattress Options	●	●
 Falls Prevention	<div><div></div></div>	<div><div></div></div>
Mobi-Lift®	●	●
Multizone Bed Exit Alarm	○	○
Secure Siderails	●	●
 Improve Patient Experience	<div><div></div></div>	<div><div></div></div>
Patient Integrated Control Panel	●	—
USB Integration	●	—
Mobi-Grips	●	●
Large Mattress Platform	●	●
Head Support	○	○

● Standard ○ Optional — Not available

SPRINT 200 E SPRINT 200

 Easy Patient Transportation	<div><div></div></div>	<div><div></div></div>
i-Drive Power® with SafetySense®	○	○
FlexiDrive®	○	○
EasyBrake®	●	●
IV&Drive®	●	●
 Simplify Caregivers Procedures	<div><div></div></div>	<div><div></div></div>
Electric Positioning (Caregiver Integrated Controls)	●	—
Scales	○	○
X-ray Cassette Holder	—	—
X-ray Sleeve (Mattress)	○	○
CPR	Manual + Electric*	Manual
ZeroGap Siderails	●	●
TR pedal at Head Side	●	●
 Hospital Operations Efficiency	<div><div></div></div>	<div><div></div></div>
Durable Parts	●	●
Removable Mattress Platform Parts	●	●
Advanced Accessories Compatibility (Head Support & Foot Support)	●	●

● Standard ○ Optional — Not available * Flattens the mattress platform


ACCESSORIES



Foot Support*



Head Support



IV Poles*



3-in-1 Shelf*



Oxygen Bottle Holder



Paper Roll Holder*

* These accessories cannot be used in conjunction with foot handles.

TECHNICAL SPECIFICATIONS

	SPRINT 200 E	SPRINT 200
Safe Working Load	320 kg / 705 lb	320 kg / 705 lb
Maximum Patient Weight	280 kg / 617 lb	280 kg / 617 lb
Average Stretcher Weight	171 kg / 377 lb	161 kg / 355 lb
Maximum Backrest / Thighrest / Calfrest Angle	90°/ 40°/ 25°	90°/ 40°/ 25°
Trendelenburg / Reverse Trendelenburg Position	+17° / -17°	+17° / -17°
External Dimensions in Standard Bed Position (length × width)	216 cm (85 in) × 89 cm (35 in)	216 cm (85 in) × 89 cm (35 in)
Mattress Platform Height	53 cm (21 in) – 86 cm (34 in)	53 cm (21 in) – 86 cm (34 in)
Castor Diameter	20 cm (7.9 in)	20 cm (7.9 in)
Clearance of Undercarriage	10.7 cm (4.2 in)	10.7 cm (4.2 in)

REFERENCES

1. LINET White paper. 2022 Best practice for patient transport in Sprint 200 for reducing WRMSDs.
2. LINET White paper. Enhanced features of Sprint 200 stretcher that support the prevention of pressure injuries. 2022



Members of LINET Group

LINET Americas, Inc.
Phone: 877-815-9895 | e-mail: info@linetamericas.com | www.linetamericas.com
©2025 Linet Americas, INC. 106092022_V

