


The Solido 3 overbed table is a practical accessory that makes the patient's stay more pleasant. It has a number of ergonomic features which make handling easy and intuitive for both the patient and caregiver.

ERGONOMIC CONTROL

Thanks to the ergonomic features and special column structure, handling of the table is effortless.

- 1 Built-in handle for easy maneuverability.
- 2 Thanks to the enhanced design, the patient can pull the table even closer to the body, supporting easier eating during mealtimes.
- 3 Exceptionally smooth adjustment of the table height, via the one-touch lever.
- 4 Low undercarriage of 1.4 in fits under any LINET bed.



SIMPLE CONTROL

The Solido 3 overbed table comes standard with a safety feature for upward direction control without the use of the control handle.

TABLETOP OPTIONS



HPL tabletop



Moulded tabletop (available in White only)

COLOR OPTIONS



White

Light Gray

Natural Maple



Beech Provence

Cherry Amber

Moka Deep Walnut

ACCESSORY



Basket

TECHNICAL PARAMETERS

Inner dimensions of eating board – Tray only (Moulded version)	14.3×13.8 in
Inner dimensions of eating board – Tray only (Moulded version)	25.2×13.8 in
Overall dimensions	37.2×17.7 in
Height adjustment	40.1×27.6 in
Load-bearing capacity	66 lbs
Undercarriage Dimensions	31.3×17.1 in
Undercarriage Height	1.4 in

LINET

LINET Americas, Inc.

Phone: 877-815-9895 | e-mail: info@linetamericas.com | www.linet.com

©2024 Linet Americas, INC. 106092022_V

LINET®, Vibwife®, Ave2® are registered trademarks of LINET Group SE



www.linet.com/en-US