

MULTICARE X

LINET

Turn to Better Care



Premium ICU Bed Designed by Nurses for Nurses

REVOLUTIONARY FRAME-BASED LATERAL TILT

CHANGING LIFE FOR THE BETTER IN INTENSIVE CARE



SAFER PATIENT HANDLING

LATERAL TILT
FOOT PEDALS

Physical strain reduction on the caregiver by 67%* mitigating risk of MSDs.

EARLY MOBILIZATION

LATERAL TILT | Mobi-Lift® | ONE-BUTTON FUNCTIONS (CARDIAC CHAIR, VASCULAR LEG)

28% shorter hospital stay* and hospital cost savings.

✓

67%
← kg →

Save 33%* of physical effort or one caregiver

✗

100%
← kg →

Standard patient transfer



EASIER PATIENT TRANSFER

LATERAL TILT

Easier patient transfer to a stretcher requiring 33%* less physical effort or one caregiver less who can attend to other duties at the same time.

PRECISE POSITIONING

1° LATERAL TILT

Gradual lateral tilt by 1° at time to gently and precisely position the patient (microshifting).

ENHANCED RESPIRATORY OUTCOMES



AUTOMATIC LATERAL THERAPY

40% functional residual capacity gain* as confirmed by Electrical Impedance Tomography.

Precise and individual adjustment of the lateral tilt based on patient's clinical needs allows to change the ventilation distribution and **helps with lung collapse.**



PREMIUM PRESSURE CARE IN THE ICU

VIRTUOSO® PRO

3-cell alternating mattress replacement with individual cell deflation.



OPTICARE® X

Fully integrated optimal immersion mattress with automatic modes.



COMBINED TILTING



Reduced risk related to repositioning patients while carrying out daily tasks with them.

SUPERIOR INTUITIVENESS & EASE OF USE



MULTIBOARD X – SIMPLIFIED AND EASY-TO-USE CONTROLS

Intuitive adjustment of the setting of both the bed and integrated mattress **from one screen**, including positioning, automatic lateral therapy setting, precision scales or patient lock controls.



FALLS PREVENTION AND SAFETY

MULTIZONE BED EXIT ALARM | i-BRAKE® | LINIS SAFETYPORT®

Lowers likelihood of patients falling by **monitoring and predicting their movements** and by **automatically breaking** the bed. Additionally giving remote overview of patients by analyzing bed usage data with **LINIS SafetyPort**.

HYGIENE CONTROL

COLUMNS CONSTRUCTION
SEAMLESS SURFACE

Best practice approach to the **prevention of hospital acquired infections**.

LINET spol. s r.o. | Želevčice 5 | 274 01 Slaný | Czech Republic
tel.: +420 312 576 400 | fax: +420 312 522 668 | e-mail: info@linet.com | www.linet.com

LINET | Art. n° 9601-1096 | Edition 08/2023 | Slight colour differences are possible. Change of technical parameters reserved.
* References for clinical claims available on the Multicare X product page at www.linet.com



Discover more